


BERLINGER SPECIAL AG



The best protection for your blood test samples

BEREG-KIT-Small

Doping control equipment for your blood sample collection on whole blood or blood serum

- Well accepted system
- Secure and reliable
- One kit - one number
- Easy to handle

BEREG-KIT-Small – security for your samples

Well accepted system

The BEREG security sealing system has been successfully used since 1994 by many anti-doping organisations, International Sports Organisations and at major events such as Olympic Games, International Championships and World Cups.

The main differences are the size (tubes up to 10ml) and the transparent plastic material. All other qualities, such as the high security level, easy handling and the clear and unique numbering system remain the same. This goes in line with the BEREG urine doping control equipment (BEREG-KIT) and rises the confidence of all involved parties.

Secure and reliable

The main security features:

- security caps
(once closed, reopening only by destroying the cap)
- the bottle is the centre of the security
- one kit – one number – corresponding numbers embossed on cap and container
- pre-sealed box with void-tape
- divers tubes
(up to 10ml) are fixed firmly inside the container (no breakage)
- tube can be centrifuged inside closed BEREG-KIT-Small



The closure system is certified by the independent Laboratory Sulzer Innotec AG, Switzerland:

«The doping control bottle BEREG-Blood-KIT provides with reasonable certitude a secure and tamperevident device that effectively prevents, protects, deters and hinders access to the bottle content.»

(date: October 12th, 2005)

Easy to handle

- only a few components needed to achieve the high security
- one container for diverse tube volumes (up to 10ml)
- can be centrifuged at site (using a table-top-centrifuge)
- compact packaging for transport and distribution
- easy to seal
- robust and resilient
- cost effective

Flexibility

One of the well-known qualities of the Berlinger Special company is their flexibility to adapt to individual customer requirements. For example any accessories can be purchased upon request



Blood-Mini-KIT

For blood screening tests on whole blood

The Blood-Mini-KIT is designed to enable an easy and secure process for screening tests on site.



Specifications of the Blood-Mini-KIT

Contents:

- 1 BD vacutainer EDTA 3ml (glass)
- 1 BD security needle
- 4 small labels (number & barcode)
- 1 security void tape (for resealing)

Dimension of kit packed in a uniquely numbered plastic box:

140 x 80 x 170mm

Weight of kit:

approx. 50g

Attention:

This kit is not designed to be sent to the laboratory – only to be used for on site screening!

Specifications of the BEREG-KIT-Small

Contents:	- 2 plastic containers with security caps - A-container with unique number in red - B-container with unique number in blue - 6-digit numbers embossed - 1 absorbent pad - packed in a styrofoam box with removable rack
Container Capacity:	maximm 10ml tubes
Dimension of styrofoam box:	140 x 80 x 170mm
Weight of kit:	approx. 140g
Accessories:	upon request

Other products:

BEREG-KIT for urine collection




BEREG-Horse-KIT for urine or blood collection




Cool box for blood sample transportation:

If you are interested in temperature controlled cooling transport systems we will be pleased to give you assistance.



Q-tag - temperature monitor

The ideal and cost efficient solution for temperature monitoring during transportation of your samples to the laboratory etc.



Developed and produced exclusively by

BERLINGER SPECIAL AG	Tel	+41 (0) 71 982 88 11
Mitteldorfstrasse 2	Fax	+41 (0) 71 982 88 39
CH-9608 Ganterschwil	E-Mail	info@berlinger.ch
Switzerland	WEB	www.berlinger.com

Berlinger Special AG - your reliable partner