

BERLINGER SPECIAL AG



The best protection for your drug test samples

BEREG-Horse-KIT

Doping control equipment for urine and blood sample collection

- Well accepted system
- Secure and reliable
- One kit - one number
- Easy to handle

BEREG-Horse-KIT – security for your samples

Well accepted system

The BEREG security sealing system has been successfully used since 1994 by many anti-doping organisations, International Sports Organisations and at major events such as Olympic Games, International Championships and World Cups.

The BEREG-Horse-KIT is a further development of the well established BEREG-KIT. It is made of plastic and has been adapted to better suit the equestrian environment.

Secure and reliable

The main differences are the size (bottles up to 175ml) and the plastic material. All other qualities, such as the high security level, easy handling and clear and unique numbering system remain the same as our human testing kit.

The main security features:

- security cap
(once closed, reopening only by destroying the cap)
- the bottle is the centre of the security
- laser engraved numbers on bottle
- corresponding number embossed on cap
- shrink-sleeved bottles
- pre-sealed box with void-tape
- firmly closed – no leakage possible



This closure system is certified by the independent Laboratory Sulzer Innotec AG, Switzerland:

«The doping control bottle BEREG-Horse-KIT provides with reasonable certitude a secure and tamper evident device that effectively prevents, protects, deters and hinders access to the bottle content.»

(date: October 12th, 2005)

One Kit – one number

The following components of the BEREG-Horse-KIT are uniquely and identically numbered:

- cap
- bottle
- box (incl. barcode)
- additional labels (incl. barcode)



Easy to handle

- can be used for either urine or blood sample collection
- compact packaging for transport and distribution
- only a few components needed to achieve the high security
- easy to seal
- robust and resilient
- cost effective



Flexibility

One of the well-known qualities of the Berlinger Special company is their flexibility to adapt to individual customer requirements.

For example any accessories can be purchased upon request. Extra kit customisation can be accommodated, e.g. extra bottles, security transport bags etc.



Specification of the BEREG-Horse-KIT

Contents of standard kit:

- 2 plastic bottles with security caps
- A-bottle with unique number in red
- B-bottle with unique number in blue
- caps with identical black number
- all numbers embossed; 7-digits
- 2 small plastic bags for transport requirements
- 1 absorbent pad
- packed in cardboard box

Bottle capacity: 175ml

Dimension of kit
packed in flat box: 185 x 115 x 84mm

Weight of kit: approx. 230g

Accessories: upon request

Other products:

BEREG-KIT for urine sample collection



BEREG-KIT-Small for blood sample collection



Developed and produced exclusively by

BERLINGER SPECIAL AG	Tel	+41 (0) 71 982 88 11
Mitteldorfstrasse 2	Fax	+41 (0) 71 982 88 39
CH-9608 Ganterschwil	E-Mail	info@berlinger.ch
Switzerland	WEB	www.berlinger.com

Berlinger Special AG – your reliable partner