

Cooling boxes

Universal cool box small / Art.-no. 94-1130

Amount of kits per box 18 BEREK-KIT small container (9 kits)
4-6 BEREK-KIT urin bottles (2-3 kits)

Performance: Summer 3°C – 9°C für 48 Stunden
Winter 2°C – 8°C für 48 Stunden
Test and performance reports are available on request.
Specific tests should be done by your own before usage.

Dimension: External :
31x21x22cm = 308x203x220mm
Internal
22x12x6.6cm = 216x120x66mm

Weight: 3.5kg per piece

Shelf life 18 months after production

Universal cool box "small" (art-no. 94-1130)
Additional cover "small" (art-no. 94-1131)
Textile outer bags "small" (separately) up on request
(10 Clip seals)



Universal cool box "large" / Art.-no. 94-1132

Amount of kits per box 30 BEREK-KIT small containers (15 kits)
6-8 BEREK-KIT urine bottles (3-4 kits)

Performance: Summer 2°C – 8°C für 48 Stunden
Winter 1.5°C – 3°C für 48 Stunden
Test and performance reports are available on request.
Specific tests should be done by your own before usage.

Dimension: External
36x31x27cm = 356x304x262mm
Internal
24x22x12.5cm = 238x213x121mm

Weight: 5.4kg per piece

Shelf life 18 months after production

Universal cool box "large" (art-no. 94-1132)
Additional cover "large" (art-no. 94-1133)
Textile outer bags "large" (separately, 10 Clip seals) up on request

